

Anti-Helicobacter pylori ELISA (IgG)

Test instruction

















ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2080-9601 G	Helicobacter pylori	IgG	Ag-coated microplate wells	96 x 01 (96)

Indication: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human antibodies of the immunoglobulin class IgG against *Helicobacter pylori* in serum or plasma for the diagnosis of infections with *Helicobacter pylori*, chronic active gastritis (antral gastritis), peptic ulcer and duodenal ulcer.

Application: The determination of specific antibodies against *Helicobacter pylori* using Anti-*Helicobacter pylori* ELISA constitutes an important and affordable starting investigation which is least stressful for the patient. Borderline or positive ELISA findings can be confirmed with the Anti-*Helicobacter pylori* EUROLINE-WB. Moreover, further information on the virulence of the pathogen can be obtained.

Principles of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with *Helicobacter pylori* antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
8. Sample buffer ready for use	light blue	1 x 100 ml	
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
 Lot description			 Storage temperature
 In vitro diagnostic medical device			 Unopened usable until



Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to +37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light ☀. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid contact with the skin.

Preparation and stability of the serum or plasma samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer.

For example: dilute 10 µl sample in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The calibrators and controls are prediluted and ready for use, do not dilute them.



Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step)

Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Mode").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (>10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values.

Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values. Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) (protect from direct sunlight).

Stopping:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open-system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1 to 4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7 to 10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. Therefore, the number of tests performed can be matched to the number of samples, minimising reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of the calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive



Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use “point-to-point” plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as “>200 RU/ml”. It is recommended that the sample be retested in a new test run at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range of non-infected persons (**cut-off value**) recommended by EUROIMMUN is **20 relative units (RU)/ml**. EUROIMMUN recommends interpreting results as follows:

<16 RU/ml:	negative
≥16 to <22 RU/ml:	borderline
≥22 RU/ml:	positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends retesting the samples.

A negative serological result does not exclude an infection. Particularly in the early phase of an infection, antibodies may not yet be present or are only present in such small quantities that they are not detectable. In case of a borderline result, a secure evaluation is not possible. If there is a clinical suspicion and a negative test result, we recommend clarification by means of other diagnostic methods and/or the serological investigation of a follow-up sample. A positive result indicates that there has been contact with the pathogen. In the determination of pathogen-specific IgM antibodies, polyclonal stimulation of the immune system or antibody persistence may affect the diagnostic relevance of positive findings. Significant titer increases (exceeding factor 2) and/or seroconversion in a follow-up sample taken after 7 to 10 days can indicate an acute infection. To investigate titer changes, sample and follow-up sample should be incubated in adjacent wells of the ELISA microplate within the same test run. For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.



Test characteristics

Calibration: As no international reference serum exists for antibodies against *Helicobacter pylori*, the calibration is performed in relative units (RU/ml).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is a bacterial lysate from the *Helicobacter pylori* strain "ATCC43504".

Linearity: The linearity of the Anti-*Helicobacter pylori* ELISA (IgG) was determined by assaying at least 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was >0.95. The Anti-*Helicobacter pylori* ELISA (IgG) is linear at least in the tested concentration range (2 RU/ml to 200 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-*Helicobacter pylori* ELISA (IgG) is 0.7 RU/ml.

Cross reactivity: The quality of the antigen used ensures a high specificity of the ELISA. Sera from patients with infections caused by various agents were investigated with the Anti-*Helicobacter pylori* ELISA (IgG).

Antibodies against	n	Anti- <i>Helicobacter pylori</i> ELISA (IgG) positive
Adenovirus	12	0%
Chlamydia pneumoniae	7	0%
CMV	11	0%
EBV-CA	12	0%
HSV-1	11	0%
Influenza virus A	12	0%
Influenza virus B	12	0%
Measles virus	12	0%
Mumps virus	12	0%
Mycoplasma pneumoniae	11	0%
Parainfluenza virus Pool	12	0%
RSV	12	0%
Rubella virus	12	0%
Toxoplasma gondii	4	0%
VZV	12	0%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.



Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 samples. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

Intra-assay variation, n = 20		
Sample	Mean value (RU/ml)	CV (%)
1	144	3.2
2	179	3.2
3	207	3.1

Inter-assay variation, n = 4 x 6		
Sample	Mean value (RU/ml)	CV (%)
1	141	3.4
2	176	3.7
3	212	4.6

Sensitivity and specificity:

59 clinically pre-characterised patient samples (INSTAND) were investigated with the EUROIMMUN Anti-Helicobacter pylori ELISA (IgG). The sensitivity amounted to 100%, with a specificity of 100%. Borderline results were not included in the calculation.

n = 59		INSTAND		
		positive	borderline	negative
EUROIMMUN Anti-Helicobacter pylori ELISA (IgG)	positive	35	2	0
	borderline	0	0	0
	negative	0	0	22

Correlation of the ELISA with a Helicobacter Urease Test: Sera from 26 patients whose bioptical material was positive in an Helicobacter Urease Test were investigated with the EUROIMMUN Anti-Helicobacter pylori ELISA (IgA and IgG).

n = 26		ELISA positive		
		IgA	IgG	IgA or IgG
Helicobacter Urease Test positive		65%	100%	100%

Reference range: The levels of the anti-Helicobacter pylori antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 20 RU/ml, 24.6% of the blood donors were anti-Helicobacter pylori positive (IgG).

Clinical significance

Helicobacter pylori (synonyms: Campylobacter pylori and C. pyloridis) is the causative agent of the second most common bacterial infectious disease in humans. It takes a mostly chronic course, with a worldwide annual death rate of 700,000 cases through stomach cancer, and affects men three times more frequently than women. Helicobacter pylori was discovered in 1983 by B. Marshall and J. Robin Warren of Perth, Western Australia. They also succeeded in cultivating Helicobacter pylori from the stomach mucous membranes of patients with chronic gastritis.

Helicobacter pylori, the only globally occurring human-pathogenic species of the genus Helicobacter, is a gram-negative bacterium of spiral form with an extremely high urease production and lives intracellularly on the luminal side of the epithelial cells of the stomach mucous membrane. Two morphological forms of the bacterium exist; a spiral form proven to be infectious, and a long-lived, coccoid form.

Humans and animals shed the bacteria in their faeces. H. pylori has the ability to survive in water. With an average worldwide prevalence of 50%, the infection rate in developing countries is considerably higher (up to 70%) than in industrial countries. In Germany, a total of about 33 million persons are infected with H. pylori, of whom about 10 to 20% develop a peptic ulcer. Antibodies to H. pylori occur in about 70% of patients with chronic active gastritis, and in 60 to 90% of cases are associated with ulcerous conditions.



Biopsies mainly show a colonisation in the mucous membrane of the antrum. In rare cases, *H. pylori* can be found in the corpus mucous membranes or, in cases of gastric metaplasia, in the duodenum. *H. pylori* infections do not heal spontaneously. The agent can persist lifelong. *H. pylori* are considered to be the causative agent of type B chronic gastritis. The clinical manifestation is dyspepsia with corresponding pain in the upper abdomen. Histologically, a superficial gastritis with atrophy of the mucous membranes might be visible. However, the large majority of infections take a clinically asymptomatic course.

With 1.65 million base pairs, the relatively small genome of *H. pylori* codes for 1,500 proteins. Some of these are special virulence factors and adhesion proteins which are characteristic for the individual *H. pylori* strains. Isolates of *Helicobacter pylori* can be broadly divided into two types, the spiral strain of type I, which can express the specific antigen protein p95 (VacA – vacuolating cytotoxin) and its associated protein p120 (CagA – cytotoxin associated antigen), as well as *H. pylori* strains of type II, which do not synthesise these proteins. Other, less toxic *H. pylori* proteins are, for example, protein p33, protein p30 (OMP, outer membrane protein) and protein p19 (OMP). As the CagA protein itself can be involved in the genesis of the tumour, it is also referred to as a bacterial oncoprotein. Alongside gastritis and the formation of ulcers, the clinical symptoms include possible late consequences, such as MALT lymphomas and adeno-carcinomas. Infections with the type I agent seem to be associated with a higher pathogenicity.

The direct detection of *H. pylori* involves taking samples from the lower third of the stomach; detection is then possible in the subsequent investigation by microscopy. The sample can also be tested for the presence of urease, and therefore indirectly for the presence of *H. pylori*. In the case of a recent infection, it is also possible to determine *H. pylori* in the faeces. *H. pylori* can be detected indirectly with a high degree of probability using a breath test, but this requires special laboratory equipment.

After contact with *Helicobacter pylori*, IgA, IgG and IgM antibodies to *Helicobacter pylori* can appear in the serum:

IgA antibody titers are mostly found after a few weeks and can still be detected over a considerable period. A positive IgA result correlates well with the activity of the gastritis. IgA antibodies are formed locally, but are not detected in the serum in every case.

IgG antibody titers frequently cannot be found until after the IgM titer has fallen. An increase in the serum IgG titer is an indication of an ongoing infection. IgG antibody titers can persist over many years. Elevated IgG antibody titers are considered to be a marker for a chronic infection.

IgM antibodies are found within a few days after contact with *Helicobacter pylori*. Specific IgM antibodies are no longer found after just a few weeks.

The persistence of *H. pylori* promotes relapses through recolonisation by residual agents from the crypts of the mucous membrane. A complete and permanent eradication of the bacteria in a diagnosed *H. pylori* infection in children, adolescents and adults leads to a reduction in the recurrence rate of 80% in the case of peptic ulcers and 20% in duodenal ulcers. The test for specific IgG antibodies against *Helicobacter pylori* is a suitable indicator for the complete eradication of the agent as a part of therapy monitoring.

Table: **Test characteristics**
of the recommended methods for the diagnosis of an *H. pylori* infection

Invasive methods (requiring gastroscopy)	Sensitivity (%)	Specificity (%)
Culture	70-90	100
Histology	80-98	90-98
Urease rapid test	90-95	90-95
PCR	90-95	90-95
Non-invasive methods (not requiring gastroscopy)		
Urease Breath Test (UBT)	85-95	85-95
Stool antigen test (using only monoclonal antibodies)	85-95	85-95
IgG antibody detection in serum	70-90	70-90



In comparison with other test methods, serological tests are less stressful for the patient. It must also be considered when performing and interpreting the various test methods that a low colonisation density (even if only transient) can lead to false-negative results in the UBT, in the stool antigen test and in the culture or histology, with the exception of serology. Low colonisation densities mainly occur under treatment with protein pump inhibitors (PPI), antibiotics, after a partial gastric resection, in cases of mucous membrane atrophy as well as in stomach cancers and MALT. The guidelines therefore recommend that for reliable *H. pylori* diagnostics, the following minimum times should be observed without *H. pylori*-suppressive treatment: two weeks after termination of PPI treatment and four weeks after a preceding eradication or other antibiotic treatment.

Hence, the serological investigation of antibodies to *Helicobacter pylori* using the Anti-*Helicobacter pylori* ELISA (IgA, IgG, IgM), the Anti-*Helicobacter pylori* IIFT (IgA, IgG or IgM) or the Anti-*Helicobacter pylori* EUROLINE-WB (IgA, IgG), by means of which cytotoxins can be clearly differentiated from other disease-relevant proteins, represents the diagnostically most reliable and inexpensive method for the diagnosis of *Helicobacter* infections and for confirming the success of treatment after eradication, and is much less stressful for the patient, particularly for children.

Literature references

1. Espinola-Klein C, Rupprecht HJ, Blankenberg S, Bickel C, Kopp H, Victor A, Hafner G, Prellwitz W, Schlumberger* W, Meyer J. (*EUROIMMUN AG). **Impact of Infectious Burden on Progression of Carotid Atherosclerosis.** Stroke 33 (2002) 2581-2586.
2. Espinola-Klein C, Rupprecht HJ, Blankenberg S, Bickel C, Kopp H, Rippin G, Victor A, Hafner G, Schlumberger* W, Meyer J; AtheroGene Investigators. (*EUROIMMUN AG). **Impact of Infectious Burden on Extent and Long-Term Prognosis of Atherosclerosis.** Circulation 105 (2002) 15-21.
3. Kist M. **S3-Leitlinie „Helicobacter pylori und gastroduodenale Ulkuskrankheit“ - eine neue Herausforderung für die mikrobiologische Diagnostik.** Mikrobiologie 20 (2010) 41-46.
4. Kanizaj FT, Katčić M, Presecki V, Gasparov S, Colić Cvrlje V, Kolarić B, Mrzljak A. **Serum antibodies positivity to 12 Helicobacter pylori virulence antigens in patients with benign or malignant gastroduodenal diseases--cross-sectional study.** Croat Med J 50 (2009) 124-132.
5. Krausse R, Garten L, Harder T, Ullmann U, Birkner T, Doniec M, Grimm H. **Clinical relevance of CagA-specific antibodies related to CagA status of Helicobacter pylori isolates using immunofluorescence test and PCR.** Infection 29 (2001) 154-158.
6. Portorreal A, Kawakami E. **Evaluation of enzyme-linked immunosorbent assay for the diagnosis of Helicobacter pylori infection in children and adolescents.** Arq Gastroenterol 39 (2002) 198-203.
7. Reynders MB, Miendje Deyi VY, Dahma H, Scheper* T, Hanke M, Decolvenaer M, Dediste A. (*EUROIMMUN AG). **Performance of individual Helicobacter pylori antigens in the immunoblot-based detection of H. pylori infection.** FEMS Immunol Med Microbiol 64 (2012) 352-363.
8. Schwarzer A, Lottspeich C, Rüssmann H, Ossiander G, Koletzko S. **Evaluation of a novel rapid one-step monoclonal chromatographic immunoassay for detection of Helicobacter pylori in stool from children.** Eur J Clin Microbiol Infect Dis 26 (2007) 475-480.
9. Shiota S, Mahachai V, Vilaichone RK, Ratanachu-Ek T, Tshering L, Uchida T, Matsunari O, Yamaoka Y. **Seroprevalence of Helicobacter pylori infection and gastric mucosal atrophy in a country with a high prevalence of gastric cancer.** J Med Microbiol (2013) Jul 5.
10. Sjö Dahl K, Lagergren J. **Epidemiological aspects of gastric adenocarcinoma: are predictive diagnostics and targeted preventive measures possible?** EPMA J 1 (2010) 461-471.
11. Storr M. **Helicobacter pylori - Update 2013** [Article in German]. Med Monatsschr Pharm 36 (2013) 172-176 quiz 177-178.
12. Veijola LI, Oksanen AM, Sipponen PI, Rautelin HI. **Association of autoimmune type atrophic corpus gastritis with Helicobacter pylori infection.** World J Gastroenterol 16 (2010) 83-88.





