

Anti-VZV ELISA (IgG)

Test instruction


















ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2650-9601 G	Varicella zoster virus (VZV)	IgG	Ag-coated microplate wells	96 x 01 (96)

Indication: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human antibodies of the IgG class against varicella zoster virus (VZV) in serum or plasma for the diagnosis of Varicella zoster virus infections.

Application: The determination of specific antibodies is the method of choice for confirmation of suspected infections (Varicella) or reactivations (Zoster) with corresponding clinical symptoms. The determination of IgG and IgM also indicates an acute infection; in the case of reactivations, specific IgA antibodies have an important diagnostic value. In both cases however, antibodies of both class IgA and IgM may be present, so that the determination of avidity for antibodies of class IgG can be of great importance. The determination of the immune status in early pregnancy and the investigation of successful immunisation is possible by means of VZV-IgG antibodies.

Principles of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with varicella zoster virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 5000 IU/l (IgG, human), ready for use	red coloured in decreasing intensity	1 x 2.0 ml	
3. Calibrator 2 500 IU/l (IgG, human), ready for use		1 x 2.0 ml	
4. Calibrator 3 100 IU/l (IgG, human), ready for use		1 x 2.0 ml	
5. Calibrator 4 10 IU/l (IgG, human), ready for use		1 x 2.0 ml	
6. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
7. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
8. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
9. Sample buffer , ready for use	light blue	1 x 100 ml	
10. Wash buffer , 10x concentrate	colourless	1 x 100 ml	
11. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
12. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
13. Test instruction	---	1 booklet	
14. Quality control certificate	---	1 protocol	
 Lot description			 Storage temperature
 In vitro diagnostic medical device			 Unopened usable until



Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light ☀. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The controls and calibrators used have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer.

Example: Add 10 µl of sample to 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.



Incubation

For **semiquantitative analysis** incubate **calibrator 3** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1 to 4** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step)

Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.

Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Wash reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values.

Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells.

Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:

Empty the wells. Wash as described above.

Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells.

Incubate for **15 minutes** at room temperature (+18°C to +25°C), protect from direct sunlight.

Stopping the reaction:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 3	P 6	P 14	P 22			C 1	P 3	P 11	P 19		
B	pos.	P 7	P 15	P 23			C 2	P 4	P 12	P 20		
C	neg.	P 8	P 16	P 24			C 3	P 5	P 13	P 21		
D	P 1	P 9	P 17				C 4	P 6	P 14	P 22		
E	P 2	P 10	P 18				pos.	P 7	P 15	P 23		
F	P 3	P 11	P 19				neg.	P 8	P 16	P 24		
G	P 4	P 12	P 20				P 1	P 9	P 17			
H	P 5	P 13	P 21				P 2	P 10	P 18			

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 4), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. Therefore, the number of tests performed can be matched to the number of samples, minimising reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 3. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 3}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive



Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 4 calibrators against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction for a patient sample lies above the value of calibrator 1 (5000 IU/l), the result should be reported as ">5000 IU/l". It is recommended that the sample be retested at a dilution of e.g. 1:400. The result in IU/l read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the reference range of non-infected persons (**cut-off value**) recommended by EUROIMMUN is **100 international Units (IU/l)**. EUROIMMUN recommends interpreting results as follows:

<80 IU/l:	negative
≥80 to <110 IU/l:	borderline
≥110 IU/l:	positive

Evaluation information: For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another EUROIMMUN recommends to retest the sample.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class IgG) can be helpful.

A negative serological result does not exclude an infection. Particularly in the early phase of an infection, antibodies may not yet be present or are only present in such small quantities that they are not detectable. In case of a borderline result, a secure evaluation is not possible. If there is a clinical suspicion and a negative test result, we recommend clarification by means of other diagnostic methods and/or the serological investigation of a follow-up sample. A positive result indicates that there has been contact with the pathogen. In the determination of pathogen-specific IgM antibodies, polyclonal stimulation of the immune system or antibody persistence may affect the diagnostic relevance of positive findings. Significant titer increases (exceeding factor 2) and/or seroconversion in a follow-up sample taken after 7-10 days can indicate an acute infection. To investigate titer changes, sample and follow-up sample should be incubated in adjacent wells of the ELISA microplate within the same test run. For diagnosis, the clinical picture of the patient always needs to be taken into account along the serological findings.



Test characteristics

Calibration: The calibration is performed in international units (IU) using the international reference serum W1044 (Anti-Varicella Zoster Serum, Central Laboratory of the Netherlands Red Cross Blood Transfusion Service, Amsterdam, Netherlands). The W1044 serum contains 50 IU/ampoule by definition and was resuspended in a concentration of 50 IU/ml.

For every group of tests performed, the extinction values of the calibrators and the international units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen used are highly purified VZV proteins (strain "Ellen") from infected, human fibroblasts (NHDF).

Linearity: The linearity of the Anti-VZV ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95 . The Anti-VZV ELISA (IgG) is linear at least in the tested concentration range (12 IU/l - 5000 IU/l).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-VZV-ELISA (IgG) is 3 IU/l.

Cross reactivity: The quality of the antigen used ensures a high specificity of the ELISA. Sera from patients with infections caused by various agents were investigated with the EUROIMMUN Anti-VZV ELISA (IgG).

Antibodies against	n	Anti-VZV ELISA (IgG)
Adenovirus	12	0%
Bordetella pertussis PT	7	0%
Bordetella FHA	12	0%
Chlamydia pneumoniae	9	0%
CMV	12	0%
EBV-CA	12	0%
Helicobacter pylori	7	0%
HSV-1	12	0%
Influenza virus A	12	0%
Influenza virus B	12	0%
Measles virus	12	0%
Mumps virus (AT)	12	0%
Mycoplasma pneumoniae	12	0%
Parainfluenza virus Pool	12	0%
Parvovirus B19	9	0%
RSV	12	0%
Rubella virus	12	0%



Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Serum	Mean value (IU/l)	CV (%)
1	269	2.3
2	335	2.5
3	2442	6.1

<i>Inter-assay variation, n = 4 x 6</i>		
Serum	Mean value (IU/l)	CV (%)
1	347	5.4
2	383	4.4
3	2983	4.9

Specificity and sensitivity: A panel of 333 clinically characterised patient samples (interlaboratory test samples from INSTAND, Labquality Finland, UK NEQUAS, MQ Schweiz) was examined with the EUROIMMUN ELISA. The test shows a sensitivity of 100% and a specificity of 100%.

n = 333		INSTAND, Labquality, NEQAS, MQ Schweiz		
		positive	borderline	negative
EUROIMMUN Anti-VZV ELISA (IgG)	positive	291	1	0
	borderline	4	1	1
	negative	0	0	35

Correlation study: 170 patient samples were investigated with the EUROIMMUN Anti-VZV ELISA (IgG) and the Wampole Anti-VZV ELISA (IgG). The EUROIMMUN ELISA showed a sensitivity of 100% at a specificity of 95% with respect to the Wampole ELISA.

n = 170		Wampole ELISA IgG		
		positive	borderline	negative
EUROIMMUN Anti-VZV ELISA (IgG)	positive	127	1	2
	borderline	0	2	0
	negative	0	0	38

Reference range: The levels of the anti-VZV antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 100 IU/l, 97.4% of the blood donors were anti-VZV positive (IgG), which reflects the known percentage of infections in adults.

Clinical significance

Varicella zoster virus (VZV), synonyme: human-pathogenic herpes virus 3 (HHV3), is the causative agent of chickenpox (varicella) after which it establishes latency and can subsequently reactivate to cause herpes zoster (shingles) [1, 2]. The virus is strictly human [3]. Chickenpox – a very contagious disease – has traditionally been regarded as a benign, inevitable disease among children (25% in 1-4-year-olds, 43% in 5-8-year-olds, 27% in 9-18-year-olds) with typical blister-like rash of the entire skin [1, 4, 5]. Now we know varicella can also be a serious infection, even in childhood, but especially in young and older adults and during pregnancy [6, 7, 8, 9].

Zoster is the endogenous recurrence of an earlier varicella infection or the result of a reinfection with existing residual immunity [2]. The average incidence of herpes zoster in Europe is 3 per 1000 people per year in the total population and more than 10 per 1000 people per year in those aged >80 years [10]. The entire virus genome is present in the latently infected ganglia. VZV is latent in multiple ganglia along the entire human neuraxis [1].



In zoster, the rash affects the spreading area of one or several sensitive nerve roots, especially T3-L3 and N. trigeminus [2, 3]. Central nervous system (CNS) complications can follow both primary infection and reactivation of VZV [8, 11]. The more serious manifestations arise when VZV invades the spinal cord or cerebral arteries after reactivation of the virus, causing diseases such as myelitis, focal vasculopathies, and encephalitis [7, 11].

Varicella also causes serious infections during pregnancy with severe consequences of maternal varicella for the infant [6, 12, 13]. At birth, maternal infection with the VZV poses a truly life-threatening risk to the newborn. The neonatal mortality rate is up to 20-30%, if the maternal VZV-infection occurs between day 4 ante partum and day 2 post partum [14]. Patients with congenital varicella syndrome (CVS) typically show clinical symptoms such as skin lesions, neurological defects, eye diseases, and/or limb hypoplasia. In rare cases, isolated manifestations in the brain or eye have been reported [15].

Antibodies against varicella zoster viruses (IgA, IgG, IgM) can be found in the serum of almost all subjects during and after a varicella infection. They can be verified by ELISA and IIFT. IgG- and IgM-antibodies against VZV are markers to confirm suspected VZV-infections, especially VZV-IgG-antibody titers during pregnancy [14, 16, 17, 18]. IgA-titers are typical for a reinfection (zoster) [19].

In addition to classic serodiagnosis of VZV, especially IgG and IgM antibodies suggestive of acute infection, measurement of VZV-IgG-avidity provides information making it possible to distinguish exactly between acute and chronic infection, as determination of avidity in other virus infections demonstrates [20]. Therefore it is of particular interest in pregnant women. Avidity describes the binding strength of specific antibody to antigen. It was found to be low in the first phase after primary infection but then to increase over time. Based on this additional information, repeated testing and unnecessary anxiety in patients can be avoided, which also is the case in the serological diagnostic of other virus infections [19, 20, 22, 23].

VZV myelitis or VZV encephalitis are diagnosed by the determination of antibodies against VZV in CSF and serum using ELISA [11]. CNS involvement results in the intrathecal synthesis of antibodies against VZV in cerebrospinal fluid (CSF). Due to the fact that specific antibodies can pass from the serum through the blood-brain barrier into the CSF by diffusion, a relative CSF/serum quotient ($CSQ_{rel.}$, synonym: antibody specificity index) [20] is determined. The quotient is calculated from the amount of specific IgG antibodies in total CSF IgG in proportion to the amount of specific IgG antibodies in total serum IgG. During conversion the CSF/serum quotient of the pathogen-specific IgG-antibody concentrations $CSQ_{path.-spec.}$ (IgG) is put into relation to the CSF/serum quotient of the total IgG concentrations CSQ_{total} (IgG) [20]. A relative CSQ result above 1.5 indicates the production of specific antibodies in the CNS and the involvement of the CNS in the disease [24, 25, 26].

Generally, a life-long immunity develops, and this is also the case after a successful protective vaccination [17, 27, 28]. A passive immunisation with specific immunoglobulins is often given to immunocompromised seronegative people, such as tumour patients and recipients of transplants, as well as seronegative pregnant women after exposure to the virus [17, 29].

Literature references

1. Dueland AN. **Latency and reactivation of varicella zoster virus infections.** Scan J Infect Dis 100 (1996) 46-50.
2. Gilden D. **Varicella zoster virus and central nervous system syndromes.** Herpes 11 (2004) 89-94.
3. Cohrs RJ, Gilden DH, Mahalingam R. **Varicella zoster virus latency, neurological disease and experimental models: an update.** Front Biosci 9 (2004) 751-762.
4. Banz K, Wagenpfeil S, Neiss A, Hammerschmidt T, Wutzler P. **The burden of varicella in Germany. Potential risks and economic impact.** Eur J Health Econ 5 (2004) 46-53.
5. Heininger U, Baer G, Bonhoeffer J, Schaad UB. **Reliability of varicella history in children and adolescents.** Swiss Med Wkly 30 (2005) 252-255.
6. Sauerbrei A, Wutzler P. **Varicella during pregnancy. 1: Epidemiology and clinical aspects.** Dtsch Med Wochenschr 17 (2004) 1983-1986.



7. Kleinschmidt-DeMasters BK, Gilden DH. **Varicella-Zoster virus infections of the nervous system: clinical and pathologic correlates.** Arch Pathol Lab Med 125 (2001) 770-780.
8. Schaller B, Bernhard P, Graber P, Steck AJ. **Cerebellar syndrome after varicella infection without virus identification in cerebrospinal fluid – an important differential ataxia diagnosis.** [Article in German] Schweiz Rundsch Med Prax 88 (1999) 1901-1907.
9. Troya Garcia J, Espinosa de Los Monteros Garde MJ, Moreno B. **Varicella pneumonia in adult population: review of 21 cases.** [Article in Spanish] Rev Clin Esp 206 (2006) 566-569.
10. Volpi A, Gross G, Hercogova J, Johnson RW. **Current management of herpes zoster: the European view.** Am J Clin Dermatol 6 (2005) 317-325.
11. Gilden DH, Bennett JL, Kleinschmidt-DeMasters BK, Song DD, Yee AS, Steiner I. **The value of cerebrospinal fluid antiviral antibody in the diagnosis of neurologic disease produced by varicella zoster virus.** J Neur Sci 159 (1998) 140-144.
12. Sauerbrei A, Wutzler P. **Varicella during pregnancy. 2. Diagnosis, prevention and therapy.** Dtsch Med Wochenschr 24 (2004) 2045-2047.
13. Martin KA, Junker AK, Thomas EE, Van Allen MI, Friedmann JM. **Occurrence of Chickenpox during Pregnancy in Woman Seropositive for Varicella-Zoster Virus.** J Inf Dis 170 (1994) 991-995.
14. Zieger W, Friese K, Weigel M, Becker KP, Melchert F. **Varicella infection at birth.** Z Geburtshilfe Perinatol 198 (1994) 134-137.
15. Sauerbrei A, Pawlak J, Luger C, Wutzler P. **Intracerebral varicella-zoster virus reactivation in congenital varicella syndrome.** Dev Med Child Neurol 45 (2003) 837-840.
16. May G. **Serologische Diagnostik rekurrerender Varizellen-Zoster-Infektionen mit Hilfe des IgA-Antikörpernachweises.** Clin Lab 42 (1996) 761-764.
17. Gross G, Schofer H, Wassilew S, Friese K, Timm A, Guthoff R, Pau HW, Malin JP, Wutzler P, Doerr HW. **Herpes zoster guideline of the German Dermatology Society (DDG).** J Clin Virol 26 (2003) 277-289. Discussion 291-293.
18. Sauerbrei A, Farber I, Brandstadt A, Schacke M, Wutzler P. **Immunofluorescence test for sensitive detection of varicella-zoster virus-specific IgG: an alternative to fluorescent antibody to membrane antigen test.** J Virol Methods 119 (2004) 25-30.
19. RKI. **Infektionskrankheiten A bis Z.** (2006).
20. EUROIMMUN AG. Stöcker W, Schlumberger W. **Alle Beiträge zum Thema Autoimmun-diagnostik.** In: Gressner A, Arndt T (Hrsg.) **Springer Lexikon Klinische Chemie. Medizinische Labordiagnostik von A-Z.** Springer Medizin Verlag, Heidelberg 1 (2007).
21. Kneitz RH, Schubert J, Tollmann F, Zens W, Hedman K, Weissbrich B. **A new method for determination of varicella-zoster virus immunoglobulin G avidity in serum and cerebrospinal fluid.** BMC Infect Dis 4 (2004) 33.
22. Levett PN, Sonnenberg* K, Sidaway F, Shead S, Niedrig M, Steinhagen* K, Horsman GB, Drebot MA. (*EUROIMMUN AG). **Use of Immunoglobulin G Avidity Assays for Differentiation of Primary from Previous Infections with West Nile Virus.** Journal of Clinical Microbiology Vol 43 (2005) 5873-5875.
23. Levett PN, Sonnenberg* K, Shead S, Sidaway F, Niedrig M, Kreher* P, Arp* B, Erdmann* S, Jahn A, Steinhagen* K, Drebot M, Horsman GB. (*EUROIMMUN AG). **IgG avidity determinations but not IgM tests are capable to differentiate recent from previous infections with West Nile Virus.** Poster zum 21. Annual clinical virology symposium and annual meeting Pan American society for clinical virology, USA, Florida, Clearwater beach (2005).
24. Reiber H, Peter JB. **Cerebrospinal fluid analysis: disease-related data patterns and evaluation programs.** J Neurol Sci 184 (2001) 101-122.
25. Reiber H, Lange P. **Virus-spezifische Antikörper in Liquor und Serum. ELISA-Analytik und Auswertung mittels Antikörper-Index und Quotientendiagramm.** Lab Med 15 (1991) 204-207.



-
26. Reiber H, Lange P. **Quantification of virus-specific antibodies in cerebrospinal fluid and serum: sensitive and specific detection of antibody synthesis in brain.** Clin Chem 37 (1991) 1153-1160.
 27. Breuer J. **Live attenuated vaccine for the prevention of varicella-zoster virus infection: does it work, is it safe and do we really need it in the UK?** J Med Microbiol 52 (2003) 1-3.
 28. Edmunds WJ, Brisson M. **The effect of vaccination on the epidemiology of varicella zoster virus.** J Infect 44 (2002) 211-9.
 29. Kangro HO, Manzoor S, Harper DR. **Antibody Avidity Following Varicella-Zoster Virus Infections.** J Med Vir 33 (1991) 100-105.



